



The Hibbot empowers your therapeutic strategies



What it is

The Hibbot is a dynamic walking aid facilitating participation in everyday activities. It enables early mobilisation with high intensity gait training in an upright, weight-bearing position.

What it does

Unlike any other conventional walkers, the Hibbot stabilizes pelvic and trunk instability with *just enough support* for the individual child. The device stimulates activity in gluteal extensor muscles and promotes excellent postural alignment and balance, a prerequisite for optimal movement function.

A game changer

The Hibbot is the first walking aid that simulates the physiotherapists' hands-on gait training, bridging a gap between therapy intervention and daily stimulation of a child's walking skills.

Indications for using the Hibbot

Children with different types of neurological disorders, and syndroms who may learn to walk either individually or using a handheld walking aid. For children with cerebral palsy they will typically have motor function level corresponding to the GMFCS levels II-III.

A good test to assess the feasibility of the Hibbot is to let the child stand and/or walk with manual support at the pelvic. If a child can stand with just this support, the Hibbot may be the ideal walking aid.

Other indications

- no structural deformities
- motivated to stand and walk



CHILDREN WITH NEUROLOGICAL DISORDERS AND SYNDROMS OFTEN HAVE:

- Abnormal muscle tone
- Substantially reduced neuromuscular activation and strength
- Reduced or affected mobility
- Poor balance
- Poor selective motor control

Learn to stand before you can walk! Through active standing the child will increase muscle strength and balance and improve its ability to walk.

How to develop walking skills

- **1. STAND** double leg support with postural control.
- **2. WEIGHT SHIFT** progression from double to single support.
- **3. TAKE STEPS** a stance phase and a swing phase of the legs.

1 // Standing activities in the Hibbot

2 // Weight shifting activities in the Hibbot



- Let the child reach for toys
- Play with a balloon; keep it in the air
- Throw a ball between the two of you
- Throw a ball into a basket
- Play bowling (hit cans or bottels)
- Kick a ball with the feets
- Draw at a table
- Cook at a play kitchen
- Play with the child's favourite toy
- Play together with friends

Start with the brake ON. To challenge the child's balance and strength switch the brake OFF.

Participation with other children is a great stimulation to develop.

- Let the child reach for a toy
- Let the child step sideways

Support and stimulate the extension of the hips with your hands when the child shifts weight from one foot to the other

Start with the compass ON. To challenge the child's balance and strength switch the compass OFF.

Sit on a low chair with wheels in front of the child. Support and stimulate with your hands extension of the hips when the child takes steps forward.



3 // Walking activities in the Hibbot



- Let the child walk with a trolley
- Let the child chase you in the hallway
- Let the child kick a ball in the hallway towards a goal drawn at the end wall
- Go shopping in the child's favourite toystore
- Explore the playground





How to adjust the Hibbot in accordance with the child's progression

- 1. Break ON/OFF
- 2. External weigths
- 3. Postural alignment
- 4. Pelvic mobility
- 5. Lateral sway Compass ON/OFF
- 6. Removable handle
- 7. Functional progress guide
- 8. Growth progress guide

1 // Brake ON and OFF

2 // External weights





With the brake on, the Hibbot will not roll backwards. Release brake temporarily please see page 17.






The Hibbot will move backwards if the child is leaning on the pelvic brace.

The child can control backward rolling by using the hip extensors to keep alignment.

Increase or decrease external weights







When the child is increasing it's strength (able to extend the hips and knees), gradually remove weights.

3 // Postural alignment

4 // Pelvic mobility in frontal and transverse plane

In sagittal plane





Anterior or posterior body inclination for optimal pelvic alignment.





Turn the blue screw on top of the pelvic plate to increase or decrease rotation in the frontal plane.







Turn the blue screw on the back of the pelvic plate to increase or decrease rotation in the horizontal plane.

Compass ON



By lateral sway the Hibbot will turn back into the centre and the child will walk straight forward.

Compass OFF



By lateral sway or sideways steps the Hibbot will not turn back into the centre. More challenging for the child.

To gradually decrease the strength of the compass or to turn it permanently off, adjust in the black box beween the wheels when the child is not in the Hibbot.



Move the handle up to temporarily release the compass.



Move the handle down to temporarily release the reverse brake.



Move the handle sideways to temporarily release both the reverse brake and the compass



Start with minimal challenge for the child by:

- Brake ON
- Weights ON
- No pelvic mobility frontal or horizontal plane
- Compass ON

As the child gains strength and stability increase challenges by:

- Brake OFF
- Gradually remove weights
- Gradually increase pelvic mobility in frontal and/or horizontal plane
- Gradually decrease compass strength or turn it permanently OFF (in the black box)

Wheel size (inch)	Pelvic Height (cm)	Axle shift (cm)	Date
	25	7	
	26	9	
	27	10	
46 "	28	12	
16 "	29	13	
	30	15	
Art.no. 500106	31	16	
	32	18	
	33	19	
	34	21	
	35 36	22 7	
	37	9	
	38	11	
	39	12	
20 "	40	14	
	41	16	
Art.no. 500107	42	17	
	43	19	
	44	21	
	45	22	
	46	7	
	47	9	
	48	11	
24 "	49	12	
24	50	14	
	51	16	
Art.no. 500108	52	17	
	53	19	
	54	21	
	55	22	
	56	7	
	57	9	
	58	10	
	59	11	
	60	12	
29 "	61	14	
Art.no. 500109	62	15	
Art.no. 500109	63	16	
	64	17	
	65	19	
	66	20	
	67	21	
	68	22	

Brace size:	XS		S		M		L	
(cm)	Art.no 500110		Art.no. 500111		Art.no. 500112		Art.no. 500113	
Pelvic-circum	38	43	43	48	48	54	54	61

Example patient:

2 year old girl of 82 cm with a pelvic height of 32 cm will have 16" wheels shifted at 18cm back from the base line

The circumference around the pelvic is 40 cm, which means she needs XS in corset.

For more information, please contact us

- **\(\)** + 47 35 50 51 20
- info.be@madeformovement.com
- www.madeformovement.com

