Lateral Trunk Supports

Spex lateral trunk supports offer adjustability, comfort and positioning. Our anatomical pad design with new full curved Natural-Lateral® surface has been tested for optimum fitting to body shape. The supporting surface of each lateral pad includes a layer of pressure relieving foam for additional protection and comfort.

CLINICAL

- Modular lateral supports provide optimum trunk control whilst improving balance, posture and comfort
- Offset correction can be applied to an asymmetric trunk to bring better trunk alignment for enhanced function

FUNCTIONAL BENEFITS


- Super comfortable pressure-reduction padding for moulding to user shape
- Removable and washable covers
- Brackets are adjustable in all planes: height, width, depth and lateral angle for ultimate trunk orientation
- Spex brings you modular creativity! Mix and match with different pad combinations
- Provide trunk control and stability with the ability to install a lateral in multiple positions on the back support
- Approx. lateral hardware weight is from 0.15 to 0.45 kg, and approx. lateral pad weight is from 0.05 to 0.3 kg

OPTIONS

- Select from Fixed, Standard Swing-Away or Axial Swing-Away Hardware
- Bracket upgrades available:
 - Cane Mount
 - DeepContour backrest bracket
 - Extension for DeepContour (commonly used with SuperShape backrest)

AXIAL® TECHNOLOGY

Unique swing-away hardware that allows the lateral pad to rotate to correctly support/align to the user's trunk needs.



HARDWARE OPTIONS



FIXED
Width and depth adjustable



STANDARD SWING-AWAY Width and depth adjustable with vertical pivot adjustment.




AXIAL® SWING-AWAY
Width, depth, and multiple angle adjustable
by it's rotational system.
Axial technology provides the correct
support structure with patented mechanically
rotational system.



CANE MOUNT

Mounts directly to wheelchair canes



Anatomical full-curve

NATURAL-LATERAL® technology
for optimum fitting to body shape

Comfort to new levels with PRESSURE-REDUCING layer


Huge pad range to provide the exact SNUG FITTING comfort that's required (refer Page 70)

FIXED LATERAL HARDWARE




Options (Add pad, shown on Pg 70)	Depth Adjustment	Code
Fixed Bracket	55 mm	1266-2101-100
Kidz Fixed Bracket	40 mm	1266-2111-100


AXIAL LATERAL HARDWARE



Options (Add pad, shown on Pg 70)	Code
Axial Hardware	1266-2301-100
Kidz Axial Hardware	1266-2311-100

SWING-AWAY LATERAL HARDWARE





STANDARD

Options (Add pad, shown on Pg 70)	Overall length	Plate Length	Code
Swing-Away, straight bracket	175 mm	120 mm	1266-2201-100
Swing-Away, kidz straight bracket	134 mm	85 mm	1266-2211-100
Swing-Away, straight extended bracket	205 mm	150 mm	1266-2202-100

OFFSET

