

# Jay Fusion™ Cushion

The JAY® Fusion™ is an adjustable cushion designed for skin protection that features JAY® Flow™ fluid or air insert options in a package that's comfortable, stable and supportive.

### COMMON CLIENT CHARACTERISTICS

- May be at high risk of skin breakdown
- May not be able to perform independent weight shifts
- May have changing postural needs over time
- May have aggressive positioning needs
- May be in a wheelchair for extended periods of time
- May experience weight loss



Fusion Rating Scale	
Side-to-Side Stability	4.5
Forward Stability	4
Pressure Management	4.5

Component	Feature	Design Intent	Clinical Value
BASE	Pre-contoured closed cell foam base	Provides lateral and forward stability that is easy to modify	Encourages orthopedic alignment, increases sitting tolerance and accommodates changing user needs
	Pre-contoured, soft foam layer with posterior pelvic wall	Provides positioning of the pelvis and thighs and reduces fluid migration	Posterior support with comfort
INSERT	Optional JAY® Flow™ fluid pad	Consistently ensures proper fluid placement by reducing fluid migration	Reduces risk of skin breakdown
	Optional air bladder available in single valve or dual valve options	Allows customization of air levels to optimize immersion and accommodate asymmetric postures	Reduces risk of skin breakdown
INNER COVER	Inner cover with Aquaquard™ zipper and anti-wicking thread	Protects foam base from moisture build up and is easy to clean	Helps with incontinence management
COVER	X-static® silver thread stretch outer cover	Reduces surface tension and resist bacterial growth	Helps prevent bacterial growth and reduces risk of skin breakdown

Technical Specifications	
Waterlow	25+
Width Range	14" to 24"
Depth Range	14" to 24"
Weight	2.4kgs 16"x18" w/fluid)
Height	4"
Weight Capacity	136kgs (14-21" widths) 227kgs (22-24" widths)
JAY Your Way Modifications Avail.	Yes
Positioning	No

		WIDTH												
		12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"
DEPTH	12"	N/A	JYW	JYW	JYW	JYW	JYW	JYW	N/A	N/A	N/A	N/A	N/A	N/A
	13"	N/A	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	N/A	N/A
	14"	N/A	JYW	•	•	•	•	•	JYW	JYW	JYW	JYW	N/A	N/A
	15"	N/A	JYW	•	•	•	•	•	•	•	JYW	JYW	JYW	JYW
	16"	N/A	JYW	•	•	•	•	•	•	•	JYW	JYW	JYW	JYW
	17"	N/A	JYW	•	•	•	•	•	•	•	•	•	JYW	JYW
	18"	N/A	JYW	•	•	•	•	•	•	•	•	•	•	JYW
	19"	N/A	JYW	JYW	•	•	•	•	•	•	•	•	•	•
	20"	N/A	JYW	JYW	•	•	•	•	•	•	•	•	•	•
	21"	N/A	JYW	JYW	•	•	•	•	•	•	•	•	•	•
	22"	N/A	N/A	N/A	JYW	•	•	•	•	•	•	•	•	•
	23"	N/A	N/A	N/A	N/A	JYW	•	•	•	•	•	•	•	•
24"	N/A	N/A	N/A	N/A	N/A	JYW	•	•	•	•	•	•	•	

• Standard size available    JYW Size available through JAY Your Way    N/A Size not available