

# Quickie Iris



## ULTIMATE TILT-IN-SPACE MANUAL WHEELCHAIR

# Client Details

Assessor Name		Client Name	
Phone:		Funding:	
Email		Trial Approval #	
Delivery Address			

# Package & Size Options

Medium Package Includes:	Long Package Includes:	
<ul style="list-style-type: none"> <li>• True-fit growth system</li> <li>• Tie downs</li> <li>• Width adjusting system</li> <li>• Medium seat frame</li> <li>• 6" x 1.5" semi pneumatic front castors</li> <li>• 16" Rear mag wheels; pneumatic</li> <li>• 80° legrest hangers with angle adjustable footplates</li> <li>• Single post height adjustable arm supports with full length pads</li> <li>• Height adjustable push handle</li> <li>• Growable seat pan</li> <li>• Anti-tips</li> <li>• Black Opal frame colour</li> <li>• Foot operated tilt</li> <li>• Foot operated brakes</li> <li>• Medium seat depth range: 16-20"</li> <li>• Spex 4 point belt</li> </ul>	<ul style="list-style-type: none"> <li>• True-fit growth system</li> <li>• Tie downs</li> <li>• Width adjusting system</li> <li>• Long seat frame</li> <li>• 6" x 1.5" semi pneumatic front castors</li> <li>• 16" Rear mag wheels; pneumatic</li> <li>• 80° legrest hangers with angle adjustable footplates</li> <li>• Single post height adjustable arm supports with full length pads</li> <li>• Height adjustable push handle</li> <li>• Growable seat pan</li> <li>• Anti-tips</li> <li>• Black Opal frame colour</li> <li>• Foot operated tilt</li> <li>• Foot operated brakes</li> <li>• Long seat depth range: 18-22"</li> <li>• Spex 4 point belt</li> </ul>	

Medium Package (Depth 16"-20")	Long Package (Depth 18"-22")
<b>Wheelchair Width 14"</b> 1567-0014-000	<b>Wheelchair Width 14"</b> 1567-0114-000
<b>Wheelchair Width 16"</b> 1567-0016-000	<b>Wheelchair Width 16"</b> 1567-0116-000
<b>Wheelchair Width 18"</b> 1567-0018-000	<b>Wheelchair Width 18"</b> 1567-0118-000
<b>Wheelchair Width 20"</b> 1567-0020-000	<b>Wheelchair Width 20"</b> 1567-0120-000

## Specifications

Seat Width (R)		
Seat Depth (S)		
Back Height (T)		
Seat to Footplate		
Seat to Floor		

## Castor Upgrades

Short Forks for Low Seat to Floor	8 x 2" Pneumatic	8 x 2" Pneumatic with Airless Insert	8 x 1.5" Semi Pneumatic
1568-8415-100	1568-8401-010	1568-8402-020	1568-8403-010

## Rear Wheel Upgrades

16" Mag with Pneumatic Airless Type	20" 5 Spoke Mag	22" 5 Spoke Mag	24" 5 Spoke Mag
1568-01-9568-116	Pneumatic 1568-9520-110	Pneumatic 1568-9522-110	Pneumatic 1568-9524-110
	Airless Insert 1568-9520-120	Airless Insert 1568-9522-120	Airless Insert 1568-9424-120
24" Spoke Wheel with Pneumatic Tyre	24" Spoke Wheel with Pneumatic Airless Tyre	BMX Tyre Upgrade	
1568-01-9524-010	1568-01-9524-020	1460-PART-055	

## Wheel Lock

Change Wheel Lock to Push to Lock	Push to Lock as additional to footlock	Push to Lock 6" Extensions	Drum Brakes 16" Pneumatic	Drum Brakes 24" Pneumatic
1568-8451-000	1568-8455-000	1579-161930	1568-8452-16P/16S	1568-8453-24P/24S

# Footrest Hangers

<b>Flared Footplate Receivers</b> (adds 2" width to legrests)	1568-01-5151-000
--	------------------


Change to 90 Degree Front Mount Hanger (when purchased with chair)	Change to 80 Degree Front Mount Hanger (when purchased with chair)	Change to 70 Degree Front Mount Hanger (when purchased with chair)	Change to 70 Degree Extension Hanger (when purchased with chair)
1568-5090-000 (pair) 	1568-5080-000 (pair) 	1568-5070-000 (pair) 	1568-5071-000 (pair) 

Front Mount Footrest Extension Tube	Elevating Leg Rests	Articulating Leg Rests
2" Pair 1568-5122-000	14-15 1568-5034-000	14-15 1568-5044-000
4" Pair 1568-5124-000	16-17 1568-5036-000	16-17 1568-5046-000
6" Pair 1568-5126-000	18-20 1568-5038-000	18-20 1568-5048-000

# Foot Plates

Heel Loops for Footplates	Footplate Padding	Spex Impact Guards
12-14" wide: 1568-5312-14	10"-11" Wheelchair: 1415-5310-300	70 Degree Hangers 1415-5360-300
15-16" wide: 1568-5315-16	12"-13" Wheelchair: 1415-5312-300	80 & 90 Degree Hangers 1415-5361-300
17-18" wide: 1568-5317-18	14"-15" Wheelchair: 1415-5314-300	Articulating Hangers 1415-5363-300
19-20" wide: 1568-5319-20	16"-17" Wheelchair: 1415-5316-300	
	18"-22" Wheelchair: 1415-5318-300	16-17" Wide
		18-20" Wide

## Contracture Footrest

Short (10.5" - 13" Leg Lengths)	Medium (12.5" - 15" Leg Lengths)	Long (14.5" - 17" Leg Lengths)	
14" 1568-5170-14S	14" 1568-5170-14M	14" 1568-5170-14L	
16" 1568-5170-16S	16" 1568-5170-16M	16" 1568-5170-16L	
18" 1568-5170-18S	18" 1568-5170-18M	18" 1568-5170-18L	
20" 1568-5170-20S	20" 1568-5170-20M	20" 1568-5170-20L	

# Foot Plates (Continued)

Jay Foot Box - Standard Height		Spex Calf Strap Available as Padded (1) or Neoprene (0) (replace the X when ordering)	
14" 1453-5014-STD		14" 1402-6x72-017	
16" 1453-5016-STD		16" 1402-6x73-017	
18" 1453-5018-STD		17" 1402-6x74-017	
20" 1453-5020-STD		18" 1402-6x75-017	
		20" 1402-6x76-017	

## Spex Calf Panels

14-15" 1402-6314-017	
16-17" 1402-6316-017	
18-19" 1402-6318-017	
20" 1402-6320-017	
Custom 1402-6329-017	

# Backrest

Dynamic Canes - (1568-2170-000)		Reclining Back Canes	
<b>Preferred Set-Up</b>		Suits sizes 14" - 20" 1568-243356-000	
Soft			
Medium			
Hard			

# Armrest

Length Adjustable Locking Flip Up Armrests	Desk Length Armrest Pad	Quickie Iris Dual Post Flip Back Armrest Upgrade
Full length pad (back cane mounted)	When purchased with chair	1568-5622-099
1568-5651-000	1568-5532-000 (each)	

# Armrest (Continued)

Spex Comfi Padded Style Armrests		Spex Armrest Attachment Upgrade
Comfi Arm Pads Size 1 Narrow 1156-5811-410	Comfi Arm Pads Size 1 Wide 1156-5821-410	Mounting to existing wheelchair armrest mounts 1156-5880-100
Comfi Arm Pads Size 2 Narrow 1156-5812-410	Comfi Arm Pads Size 2 Wide 1156-5822-410	
Comfi Arm Pads Size 3 Narrow 1156-5813-410	Comfi Arm Pads Size 3 Wide 1156-5823-410	

# Tray

Spex Clear Tray 5ml	Spex Black Tray 5ml	Zippie Style Tray	Travel Tray
Small 1311-7111-000 (suits 10-13" wheelchairs)	Small 1311-7101-000	X-Small 1560-7110-100	Small 1311-7441-209
Medium 1311-7112-000 (suits 14-17" wheelchairs)	Medium 1311-7102-000	Small 1560-7113-100	Medium 1311-7442-209
Large 1311-7113-000 (suits 18-21" wheelchairs)	Large 1311-7103-000	Medium 1560-7115-100	Large 1311-7443-209
		Large 1560-7118-100	
		X-Large 1560-7121-100	
			

# Tray Accessories

Elbow Protraction Guards		Elbow Pads			
Guards for Small Tray 1311-7351-000		Neoprene Pads 8mm 1311-731X-400 <small>Replace the X with the code for the size required</small>	Small 1	Medium 2	Large 3
Guards for Medium Tray 1311-7352-000					
Guards for Large Tray 1311-7353-000					

## Tray Padding

	8mm Neoprene Padding	10mm CCF Foam Padding	
<b>Small Tray</b>	1311-7431-291	1311-7431-360	
<b>Medium Tray</b>	1311-7432-291	1311-7432-360	
<b>Large Tray</b>	1311-7333-291	1311-7343-360	

# Canopy Options

Iris Canopy	Folding Sunhood	UV Shade	Rain Cover
Suits 10" - 14" Chairs 1568-125320	450mm wide 1446-7545-000	1446-7710-000	2880-7651-000 (bingo buggy style)
Suits 15" - 20" Chairs 1568-125321	550mm wide 1446-7555-000		
			

## Notes

### Further Information

---

---

---

---

---

---



# medifab

Distributed by Medifab Ltd

32 Detroit Drive  
Rolleston 7675  
New Zealand

+64 3 307 9790

[solutions@medifab.com](mailto:solutions@medifab.com)

[www.medifab.com](http://www.medifab.com)



DISTRIBUTED BY:

Script: Quickietris\_20220727