





1



**MEASURE WIDTH ACROSS BACK SURFACE**  
SIMPLY MEASURE ACROSS THE BACK BETWEEN AXILLAE

USER BACK WIDTH	ZYGO WIDTH
8"-10"	10"
10"-12"	12"
12"-14"	14"
14"-16"	16"
16"-18"	18"
18"-20"	20"


2



**MEASURE FOR BACK SUPPORT HEIGHT**

SIMPLY MEASURE FROM THE DESIRED  
FINISH HEIGHT TO 1" BELOW PSIS

3



**CHOOSE CONTOUR**

**ACTIVE**




**DEEP**

ACTIVE LO

DEEP LO





ACTIVE MID

DEEP MID






**SIZING TABLE - ACTIVE LO**

ZYGO SIZE	OUTER WIDTH	HEIGHT	DEPTH	LATERAL HEIGHT
12" X 7"	10.7"	7.8"	2"	7"
14" X 8"	12.7"	8.6"	2.3"	7.5"
16" X 9"	14.7"	9.8"	2.6"	8.4"
18" X 9"	16.7"	10.3"	2.8"	8.9"
20" X 11"	18.7"	11.5"	3.3"	9.9"







**SIZING TABLE - ACTIVE MID**

ZYGO SIZE	OUTER WIDTH	INNER WIDTH	HEIGHT	DEPTH	LATERAL HEIGHT	LATERAL POSITION
10" X 9"	8.8"	N/A	9"	1.4"	4.3"	2.3"
10" X 11"	8.8"	N/A	11"	1.4"	4"	2.8"
12" X 10"	10.7"	N/A	10"	1.8"	4.8"	2.5"
12" X 12"	10.7"	N/A	11.5"	1.8"	4.9"	3.1"
12" X 13"	10.7"	N/A	12.5"	1.8"	5.3"	3.1"
14" X 13"	13.5"	N/A	13"	2.2"	6.6"	2.9"
14" X 15"	13.5"	N/A	15"	2.2"	6.1"	3.6"
14" X 16"	13.5"	N/A	16"	2.2"	6.1"	3.8"
16" X 14"	15.3"	N/A	14"	2.6"	7"	3.1"
16" X 16"	15.3"	N/A	16"	2.6"	6.3"	3.6"
16" X 17"	15.3"	N/A	17"	2.6"	7.2"	3.4"
18" X 15"	17.2"	N/A	14.5"	2.9"	7"	3.3"
18" X 17"	17.2"	N/A	16.5"	2.9"	6.6"	3.6"
18" X 18"	17.2"	N/A	18"	2.9"	6.6"	4.4"
20" X 15"	19"	N/A	15"	3.3"	7.4"	3.5"
20" X 18"	19"	N/A	17.5"	3.3"	6.7"	4.2"



**SIZING TABLE - DEEP LO**

ZYGO SIZE	OUTER WIDTH	INNER WIDTH	HEIGHT	DEPTH	LATERAL HEIGHT
12" X 7"	10.5"	8.5"	7"	3.5"	6.2"
14" X 8"	12.5"	10.5"	8.2"	4.5"	7.4"
16" X 9"	14.7"	12.5"	8.8"	5.8"	7.8"
18" X 9"	16.7"	14.5"	9.4"	6.3"	8"
20" X 11"	18.7"	16.5"	10.6"	6.3"	8.5"



**SIZING TABLE - DEEP MID**

ZYGO SIZE	OUTER WIDTH	INNER WIDTH	HEIGHT	DEPTH	LATERAL HEIGHT	LATERAL POSITION
10" X 9"	8.8"	6.8"	9"	2.5"	4.3"	2.3"
10" X 11"	8.8"	6.8"	11"	2.5"	4"	2.8"
12" X 10"	10.8"	8.8"	10"	3"	4.8"	2.5"
12" X 12"	10.8"	8.8"	11.5"	3"	4.9"	3.1"
12" X 13"	11"	9"	12.5"	3.8"	5.3"	3.1"
14" X 13"	12.9"	10.5"	13"	3.6"	6.6"	2.9"
14" X 15"	12.9"	10.5"	15"	3.6"	6.1"	3.6"
14" X 16"	13.5"	11.2"	16"	4.7"	6.1"	3.8"
16" X 14"	15.1"	12.8"	14"	5.2"	7"	3.1"
16" X 16"	15.1"	12.8"	16"	5.2"	6.3"	3.6"
16" X 17"	15.2"	12.9"	17"	6"	7.2"	3.4"
18" X 15"	17"	14.7"	14.5"	5.7"	7"	3.3"
18" X 17"	17"	14.7"	16.5"	5.7"	6.6"	3.6"
18" X 18"	17.3"	14.9"	18"	6"	6.6"	4.4"
20" X 15"	19"	16.7"	15"	6.3"	7.4"	3.5"
20" X 18"	19"	16.7"	17.5"	6.3"	6.7"	4.2"