


## Wheelchair Back Supports


### VIGOUR HI

Solid full height planar back support for maximum trunk and shoulder support, suitable for users who require minimal to moderate postural support.




### ZYGO ACTIVE

Lo and mid height shell options with an open, shallow contour and inbuilt laterals, suitable for users with moderate trunk mobility and strength.




### ZYGO DEEP

Lo and mid height shell options with a supportive, deep contour and inbuilt laterals, suitable for users with limited trunk mobility and strength.




### SPEX STRAIGHT

Contouring cube system allows on-the-spot customisations for torso comfort and asymmetry needs whilst providing greater stability for the upper body.




### SPEX T-SHAPE


T-shape, with contouring cube system to allow on-the-spot customisations for torso comfort and asymmetry needs whilst providing greater stability for the upper body.




### BACK HARDWARE




**Quad Mount**  
1259-2400-100  
4 sturdy points of attachment, default for Vigour and Spex



**Drop Mount**  
1259-1400-00 [P: Paediatric, S: Small, M: Medium, L: Large]  
Default for Zygo



**Dual Mount**  
1259-2200-100  
For increased adjustability, angle & growth



**Heavy Duty Quad Mount**  
1259-2405-100  
Default for back supports over 20" wide or 23" high

#### ADD A SPLASH OF COLOUR

Upgrade your back support with a splash of colour, for no extra cost! (unavailable for trial). We offer 9 colour options for our back support outer covers. Simply add the upgrade codes to your back support order.



- Granite Grey -090
- Ocean Blue -060
- Royal Purple -077
- Calypso Pink -028
- Chilli Red -022
- Mango Orange -030
- Bumblebee Yellow -044
- Leaf Green -054

#### OUR BACK SUPPORT CODES ARE MEANINGFUL!

Here's how to generate a product code to place your order.

**EXAMPLE**  
A Vigour Hi Back Support of width 18" and height 20", has the code 1252-1820-00H

**BACK SUPPORT STYLE**  
1 = Zygo  
2 = Vigour  
4 = Manta  
5 = Spex

**BACK SUPPORT WIDTH**  
e.g 16"

**BACK SUPPORT HEIGHT**  
e.g 18"

**BACK SUPPORT STYLE**  
AL= ACTIVE LO ZYGO  
AM= ACTIVE MID ZYGO  
DL= DEEP LO ZYGO  
DM= DEEP MID ZYGO  
OH= VIGOUR HI  
SS = SPEX STRAIGHT  
TS = SPEX T-SHAPE



- Blue squares are sizes which are available on the MOH list, available for straight purchase.
- Red squares are standard sizes available for trial.
- Non-standard or custom sizes are available for purchase.