

# Accessible list

## Rigid Shell Back Supports

### VIGOUR BACK SUPPORT

A tough & lightweight back support meeting the daily lifestyle demands of active wheelchair users.

Back support angle and seat depth offers exceptional adjustability and can be facilitated with the user in the wheelchair.

Ventilated spacer-fabric and shell system increases air-flow through the back offering daily comfort.

In combination with lateral supports offers exceptional trunk and thoracic stability.

Can be upgraded with **DIFFERENT FOAM PADDING OPTIONS** for excellent results

Dual-layer cover reduces **SHEAR FORCES** with air-mesh side panels to counter moisture and temperature build-up

Add laterals as required from our **LARGE MODULAR RANGE**

**SLOTS & HOLE POSITIONS** on shell allow easy installation of head supports, lateral supports and harnesses




### VIGOUR MID

Lightweight and designed for efficient propelling, with scapular subtractions for unimpeded arm function.

Designed for the individual who requires:


- A rigid shell back support with scapula cut-outs for self-propel users.
- Minimal to moderate postural support
- Improved spinal and pelvic positioning whilst maintaining propelling efficiency
- For clients that have mild torso asymmetries




#### SELECT YOUR HARDWARE

A range of hardware is offered for the Vigour back supports range. Quad mount is the default hardware.


Quad Mount  
1259-2400-100  
(Default)




Single Mount  
1259-2100-100  
(For back supports 18" high or lower)



Dual Mount  
1259-2200-100



Heavy Duty Quad Mount  
1259-2405-100




### VIGOUR HI

Solid full height back support available in three heights for maximum trunk and shoulder support.

Designed for the individual who requires:

- A solid back support
- Minimal to moderate postural support
- Improved spinal and pelvic positioning whilst maintaining propelling efficiency

Supports effective use of upper body harness accessories.



Blue squares are sizes which are available on the MOH list, available for straight purchase.  
Red squares are standard sizes available for trial.  
Non-standard or custom sizes are available for purchase.

#### MICRO CUBE FOAM UPGRADE

Foam upgrades allow air-flow breathability for maximum control of temperature and moisture. Micro-cube cells provide zoned weight distribution enhancing comfort.



#### NEW COLOURS

Upgrade your back support with a splash of colour, for no extra cost! (unavailable for trial). We offer 9 colour options for our back support outer covers.


Simply add the upgrade codes to your back support order.

- Leaf Green 1209-2491-054
- Granite Grey 1209-2491-090
- Royal Purple 1209-2491-077
- Calypso Pink 1209-2491-028
- Bumblebee Yellow 1209-2491-044
- Mango Orange 1209-2491-030
- Chilli Red 1209-2491-022
- Ocean Blue 1209-2491-060

#### OUR BACK SUPPORT CODES ARE MEANINGFUL!

Here's how to generate a product code to place your order.

**EXAMPLE**  
A Vigour Hi Back Support of width 18" and height 20", has the code 1252-1820-00H



#### PRODUCT DETAILS

Width	Capacity	Product Weight
10"-15"	200lbs (90kg)	Mid Back (16x16) 1.0kg
16"-20"	300lbs (136kg)	Hi Back (16x16) 1.3kg

Note: This is the weight of the back support only and does not include the hardware. For hardware weight, refer to the Spex Catalogue.

#### PRODUCT PERFORMANCE

Postural Control	4
Lateral Stability*	3
Postural Asymmetries	3
Changing Needs	3
Pressure Care	3

Note: This is a guide only and ratings are dependent on correct use of product.  
\*without laterals installed