

## QUICKIE 2



**Package includes:**

- DAB depth & angle adjustable backrest
- 24" lite spoke rear wheels
- Anodised push rims
- Marathon plus tyres
- 5" x 1.5" semi-pneumatic castors
- Single post height adjustable armrests with full length pads
- 80 degree footrest hangers
- Angle adjustable footplates
- 16" high 3DX tension adjustable back upholstery
- Integrated push handles
- Anti-tips
- Spex 2 point belt
- Colour is Black Opal

Stocked Sizes Available (Depth)				
	16"	17"	18"	20"
Width	14" 1571-1416-001			
16"	1571-1616-001		1571-1618-001	1571-1620-001
17"		1571-1717-001		
18"	1571-1816-001		1571-1818-001	1571-1820-001
20"			1571-2018-001	1571-2020-001

Seat Width	12-22"
Seat Depth	112-20"
Front Seat Height	13.5 to 22"
Rear Seat Height	13 to 22"
Frame Style	Fixed or SA front
Back Height	9-21"
Back Support Angle	-15° - +3° (Optional)
Frame / Hanger Angle	60°/ 70°/ 80°/ 90°
Inset	N/A
Camber	(Adjustable) Set at 3°
Max. User Weight	136 kg Std / 158 kg HD
Product Weight	11.5 kg
Transit Approved	Yes
Frame Colour	24 colours
Warranty	Limited Lifetime
ARTG	100167
Folding	Yes

## ZIPPIE IRIS



**Package includes:**

- True-fit growth system
- Tie downs
- Width adjusting system
- Short seat frame
- 5" x 1.5" semi pneumatic front castors
- 16" Rear mag wheels; pneumatic
- 90° legrest hangers with angle adjustable footplates
- Offset hanger mounts
- Single post height adjustable arm supports with full length pads
- Height adjustable push handle
- Growable seat pan
- Anti-tips
- Black Opal frame colour
- Foot operated tilt
- Foot operated brakes
- Spex 4 point belt

Stocked Sizes Available (Depth)			
	10-14"	12-16"	14-18"
Width	10" 1562-0010-090	1562-0110-090	
12"		1562-0112-090	1562-0212-090
14"			1562-0214-090

Seat Width	10" to 18"
Seat Depth	10" to 18"
Front Seat Height	12.5" to 19.5"
Rear Seat Height	11" to 17.5"
Frame Style	Rotation in Space
Back Height	15-24"
Back Support Angle	-5° to +30°
Frame / Hanger Angle	60°/ 70°/ 80°/ 90°/ Elevating/Z-finity
Inset	N/A
Camber	0°
Max. User Weight	102 kg
Product Weight	10 kg
Transit Approved	Yes
Frame Colour	27 colours
Warranty	Limited Lifetime on frame. 1 year on parts.
ARTG	100167
Folding	Yes

## QUICKIE IRIS



**Package includes:**

- True-fit growth system
- Tie downs
- Width adjusting system
- Medium seat frame
- 6" x 1.5" semi pneumatic front castors
- 16" Rear mag wheels; pneumatic
- 80° legrest hangers with angle adjustable footplates
- Single post height adjustable arm supports with full length pads
- Height adjustable push handle
- Growable seat pan
- Anti-tips
- Black Opal frame colour
- Foot operated tilt
- Foot operated brakes
- Spex 4 point belt

Stocked Sizes Available (Depth)		
	16-20"	18-22"
Width	14" 1567-0014-000	1567-0114-000
16"	1567-0016-000	1567-0116-000
18"	1567-0018-000	1567-0118-000
20"	1567-0020-000	1567-0120-000

Seat Width	14-22"
Seat Depth	14-22"
Front Seat Height	12.5-19"
Rear Seat Height	12.5-19"
Frame Style	Rotation In Space
Back Height	15-24"
Back Support Angle	-5° to +30°
Frame / Hanger Angle	60°/ 70°/ 80°/ 90°/ Elevating
Inset	N/A
Camber	0°
Max. User Weight	136 kg / 181 kg HD
Product Weight	17.6 kg
Transit Approved	Yes
Frame Colour	25 colours
Warranty	Limited Lifetime
ARTG	100167
Folding	Yes



Configure with Jay Seating



Configure with Spex Seating



Full range of Sunrise wheelchairs available, including:



QXI



M6



KRYPTON R



7<sup>R</sup>



ZIPPIE<sup>®</sup> XCAPE